



Note 1:

Mounting Instruction - see 504100.
 Instructions de montage - voir 504100.
 Montageanleitung - siehe 504100.
 Monteringsvejledning - se 504100.

Item	Qty	Drawing no	Description3	Rev no
1	1	502934	Console for Dahl Multi Adaptation Set	01
2	1	503085	Locking Plate Light with Spacer Complete	01
3	2	503081	Locking Socket Complete	00
4	1	503084	Cross Tube with Locks	03
5	1	503090	Cross Tube	02
6	1	502021	Adaptation Plate Complete, Right and Left	00
7	1	502022	Hook Complete, Right and Left	00
8	4	503270	Clamp with Insert Ø20 mm	00
9	1	503653	Sticker "Millimeter"	00
10	2	504069	Sticker "MADS"	00

✓	Dimension	Material:	List of spare parts:	Weight: 0 kg
	Chamfer:	Surface treatment:	Measures without tol.: DS/EN 22768-M	
Latest correction:			Replaced by:	Replaces:
		Date:	24/01/2025	
		Issuer:	LJ	
Title: Docking adaptionssæt, aftagelig version Kit d'adaptation du système d'arrimage, version détachable 48 cm Docking-Adaptionskit, abnehmbare Version 48 cm Docking Adaptation Kit, Detachable Version 48 cm		Proactiv Traveler 4all Ergo		Scale:
				Drawing no.:
				503990
				00